EXACTA

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- **Capacity of 17 litres**
- **Fully microprocessed** management for a constant monitoring of parameters and a precise control of all the functions of the sterilizer
- **Initial Vacuum through pump** that enables the efficient extraction of air allowing for optimal steam saturation and penetration even in handpieces and turbines
- **Combined vacuum drying system** (Mocom patent) that sums up the action of a vacuum suction pump and forced air ventilation pushed through a bacterial filter
- **Fast cycle** that enables sterilisation in only 14 minutes avoiding both mechanical vacuum and vacuum drying phases
- **9 pre-set cycles** plus 1 programmable and 1 special cycle
- **Dry-off push-button** that excludes the vacuum drying phase reducing cycle time of 12 minutes
- **Easy to use**
- **Separated tanks to prevent re-cycling** for constant use of pure water and to increase the life of the sterilizer
- **It can be built-into** furniture
- **Mechanical safety** and electronic devices

For each requirement, the most suitable cycle

- **1 fast cycle:**
134°C - 4 minutes (total time 14 minutes)
- **9 pre-set cycles:**
- 6 instrument sterilisation cycles**
 - 2 cycles for solid loads
134°C - 4 minutes (total time 35 minutes)
121°C - 20 minutes (total time 40 minutes)
 - 2 cycles for wrapped loads
134°C - 12 minutes (total time 43 minutes)
121°C - 30 minutes (total time 55 minutes)
 - 2 cycles with slow decompression phase
134°C - 12 minutes (total time 55/60 minutes)
121°C - 30 minutes (total time 60/65 minutes)
- 3 disinfection cycles**
 - 1 cycle for solid loads
105°C - 20 minutes (total time 41 minutes)
 - 1 cycle for wrapped loads
105°C - 30 minutes (total time 51 minutes)
 - 1 cycle with slow decompression phase
105°C - 30 minutes (total time 60 minutes)
- **Dry-off push-button:**
Reduces cycle time of 12 minutes
- **1 programmable cycle:**
Enables the setting of temperature and exposure time allowing for the sterilisation of heat-sensitive materials at lower temperatures (min. 115°C).
- **1 Special cycle: Vacuum Test**

EXACTA

Steam sterilizer

THE SAFE, FAST AND COMPETITIVE STERILIZER

- The water steam sterilizer with initial vacuum, vacuum drying and fast cycle
- Efficient air removal
- Perfect drying
- Fast cycle
- Sterilisation of thermolabile material (programmable cycle)
- Preservation of instruments against oxidation
- Provides saturated steam during the entire cycle

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
Mocom: solid experience for higher quality

Mocom is ISO 9002-EN 46002 certified and is member of the European Standards Commission (CEN - TC 102 wg 5) for small steam sterilizer

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Microprocessor, temperature and pressure sensors

Temperature and pressure sensors, capable of communicating with the microprocessor, ensure optimal conditions for the safest instrument sterilisation. In fact, the analogical sensors used by Mocom provide higher performance levels in terms of precision and stability allowing for constant crossed controls between pressure and temperature for a constant balanced saturated steam.

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Facility for printer

The Exacta sterilizer is equipped with facilities for connections to a printer, allowing minute by minute and step by step hard-copy recordings of time, temperature and pressure.

Cyclus and DPU Printers

Cyclus and DPU are the optional compact and ergonomic printers that Mocom recommends to its Customers.

Standard Kit

- 3 trays for instruments
- 1 stainless steel support
- 1 tray remover
- 1 spare bacteriological filter
- 1 water filling tank complete of 1 pipe with quick connector
- 1 waste pipe with quick connection

Optional Accessories

- Cyclus and DPU printers
- Spiral envelope holder frame

Klave simulator

A software programme developed by Mocom for the diagnosis of anomalies on PC (for technical service).

Quality System certified by:

Via Delle Azalee, 1 - 20090 Buccinasco (MI) Italy - Tel. (+39) 024570.1505 r.a. - Fax (+39) 024884.0514
www.mocom.it - e-mail: mocomcom@interzona.com