

mocom Steriliser

Cycle Times

Mocom Millennium B+ & B²

Material	Cycle choice	Total cycle time inc. drying		Cycle temp./ holding time
		B+	B ²	
Wrapped hollow* instruments, porous materials	Type B	30 mins	40 mins	134°C 4 mins
Wrapped hollow* instruments, porous materials	Type B	44 mins	54 mins	134°C 18 mins
Unwrapped hollow* instruments	Hollow	24 mins	36 mins	134°C 4 mins
Wrapped solid instruments	Wrapped	23 mins	30 mins	134°C 4 mins
Unwrapped solid instruments	Solid	17 mins	24 mins	134°C 4 mins
Unwrapped solid instruments	Emergency	12 mins	14 mins	134°C 3 mins
Wrapped hollow* instruments, porous materials	Type B	45 mins	54 mins	121°C 20 mins
Unwrapped hollow* instruments	Hollow	38 mins	50 mins	121°C 20 mins
Wrapped solid instruments	Wrapped	37 mins	44 mins	121°C 20 mins
Unwrapped solid instruments	Solid	31 mins	38 mins	121°C 20 mins

* Hollow instruments type "A" and "B" according to EN 13060:2004

Mocom Exacta S & S²

Material	Cycle choice	Total cycle time inc. drying		Cycle temp./ holding time
		Exacta S	Exacta S ²	
Wrapped solid and hollow** instruments	Wrapped	45 mins	38 mins	134°C 4 mins
Wrapped solid and hollow** instruments	Prion	59 mins	52 mins	134°C 18 mins
Unwrapped solid and hollow** instruments	Solid	34 mins	28 mins	134°C 4 mins
Unwrapped solid instruments	Emergency	26 mins	18 mins	134°C 3 mins
Wrapped solid and hollow** instruments	Wrapped	59 mins	50 mins	121°C 20 mins
Unwrapped solid and hollow** instruments	Solid	48 mins	40 mins	121°C 20 mins

** Hollow instruments type "B" according to EN 13060:2004